



INTRODUCING MIST 360°

**SAME ACCELERATED HEALING OF MIST® THERAPY
IN A NEW APPLICATOR DESIGN**

- ▶ Improve access to difficult to treat locations
- ▶ Enhance wound visualization during treatment
- ▶ Minimize repositioning of your patients
- ▶ Treat larger wounds without stopping to refill saline
- ▶ Lighter transducer designed to minimize user fatigue

TWO APPLICATOR SOLUTIONS WITH THE SAME PROVEN RESULTS



Choose the easy to use applicator that best addresses your wound challenge.

FEATURE/ FUNCTIONALITY	MIST® 360°	MIST®
		
Saline Delivery	Provide your own saline (IV saline) <ul style="list-style-type: none"> • Offers flexibility to choose the right amount of saline for the wound size • No refilling required 	100 cc bottle of saline included in kit <ul style="list-style-type: none"> • Each bottle delivers saline for approximately 5 minutes before refilling is required
Wound Access	Flexibility to treat difficult to reach locations/odd angles <ul style="list-style-type: none"> • Minimizes the need to reposition your patient 	Provides 30° of rotation <ul style="list-style-type: none"> • Great option for easy to access, smaller wounds
Packaging	Case of 12 <ul style="list-style-type: none"> • Eco-friendly pouch offers less packaging and is easy to dispose 	Case of 6 <ul style="list-style-type: none"> • Classic tray configuration functions as waste receptacle during treatment

ORDERING INFORMATION

MIST Therapy System	Each	CP-80004	
MIST Therapy Cart	Each	CP-80008	
Applicator Kit - MIST	Box of 6	CP-80011	Each kit includes applicator, 2 germicidal wipes, absorbent pad, 100ml bottle of saline
- MIST 360°	Box of 12	CP-80024	Each kit includes applicator, 2 germicidal wipes, absorbent pad*

*Other Materials Required

0.9% sodium chloride injection USP (normal saline) 100-500cc bag

500 pressure infusor for an I.V bag (gauge that measures to 300mm Hg)